

Resuscitation Trolleys

ARTEMIS 

All Artemis Medical Resuscitation Trolleys are manufactured in the UK to our own specifications.

We offer two main types of trolley: a heavy duty, custom-built unit with all metal drawers and a lighter Vari-Cart based modular unit.

Constructed of steel and nylon-coated, they are built to withstand the rigours of Intensive Care life for many years.

A range of standard designs are available to suit most circumstances. However, we do offer a customising facility should you require it.



Heavy Duty Trolley with left-hand Melamine drop-flap and rail assembly



Custom Cart



Custom Cart

Artemis Medical Ltd
Butterly Avenue, Questor Business Park,
Dartford, Kent DA1 1JG
Tel: 01322 628877 • Fax: 01322 628878
www.artemismedical.co.uk

Heavy Duty Trolleys

FW2100	Emergency Cart 800 (w) x 550 (d) x 1000 (h) mm White nylon coated, complete with: Frame including handle Top work surface 720(w) x 550(d) mm c/w cut-out for cylinder 4 x high quality 4" castors (2 with brakes), 4 x buffers Cabinet to fit into frame FW2101 Small drawer FW2102 Medium drawer FW2103 Large drawer FW2104 Extra large drawer
FW2105	Divider to fit 2101
FW2106	Divider to fit 2102
FW2107	Divider to fit 2103
FW2008	Notes box/Clipboard holder - A4 size
FW2009	Medical rail fitted to rear & side of cart
FW2010	Drip pole (right or left hand side)
FW2011	Rail assembly with drip pole (Two rails - height adjustable)
FW2012	Drop flap
FW2013/S	Single cylinder holder - E size
FW2013/D	Double cylinder holder - E size
FW2014	Container storage bracket - rail mounted Other options available on request



FW2100 with drop-flap

Rail Trolleys

FW2700	Freeway Rail Trolley 630(w) x 730(d) x 1480(h) mm - complete with 8 lengths of rail
FW2701	Rail Trolley Melamine Shelf 540 x 370mm
FW2702	Rail Trolley Shelf with cut out 540 x 690mm
FW2703	Rail Trolley Bottom Shelf 540 x 360mm
FW2704L	Rail Trolley Shelf Drop Flap 370 x 350mm, left hand
FW2704R	Rail Trolley Shelf Drop Flap 370 x 350mm, right hand
FW2705	Rail Trolley Container Storage Panel 610 x 150mm
FW2706	Green Plastic Container - small
FW2707	Green Plastic Container - large
FW2708	White Metal Drawer Unit Small - "Easy run"
FW2709	White Metal Drawer Unit Medium - "Easy run"
FW2710	White Metal Drawer Unit Large - "Easy run"
FW2710A	White Metal Drawer complete with Divider Large - "Easy run"
FW2711	Drip Stand for Rail Trolley
FW2712	Drawer Unit Mounting Bracket
FW2713L	Metal Drop Flap for Rail Trolley Left hand
FW2713R	Metal Drop Flap for Rail Trolley Right hand
FW2714	Rail Trolley Top Metal Shelf
FW2715	Rail Trolley Modular Storage Drawer Assembly complete with: 4 x FW3001 Clear Plastic Tray 1 x FW3013 Deep White Plastic Coated Basket



FW 2700
Rail Trolley with drop flap,
drip pole and modular storage
drawer assembly

Vari-Cart Trolleys

FW3100	Resuscitation Vari-Cart - single column White nylon coated cart with flush easy-clean top & 4 high quality castors (2 with brakes) Drawer runners complete with: FW3013 Basket - deep x 1 FW3001 Plastic tray - clear x 4
FW3100R	Resuscitation Vari-Cart - single column As FW3100 but red nylon coated cart
FW3001	Clear trays - high impact polystyrene. 100mm deep. Complete with integral drawer pull
FW3002	Drop flap - Melamine
FW3003	Medical rail
FW3004P	Plastic coated drawer dividers
FW3004S	Stainless steel drawer dividers
FW3005	Stainless steel top
FW3006	Anti-static castors (x 4)
FW3007	Metal drawer
FW3008	Metal bottom shelf
FW3009	Drip pole
FW3010	Notes box - A4 size
FW3011	Clipboard & holder - metal, A4 size
FW3012	Basket - shallow
FW3013	Basket - deep
FW3014	Flush Melamine top
FW3017	Drop-in top to suit FW3100
FW3020	Link brackets
FW3021	Drop flap - metal - plastic coated
FW3022	Metal tray - various colours nylon coated
FW3023	Metal tray - deep - various colours nylon coated
FW3024	Metal shelf for mounting on runners - various colours
FW3025	Vinyl cover - for FW3100