


# Transmitters

ZS-910/920/930/940PK  
ZM-940PK

**Better care  
by wireless**



Emergency

Transport

Ward

Nurse Station

*Fighting Disease with Electronics*

**NIHON KOHDEN**

# High quality patient care by wireless

—patient is free from the bed—

## ❖ Wireless monitoring system

You can make a wireless monitoring system with transmitters and a Nihon Kohden central monitor with multiple patient receiver or a Nihon Kohden telemetry system.

## ❖ Increase patient safety

By continuous measurement of important vital parameters, patient condition is always monitored so changes can be quickly handled.

## ❖ Increased patient comfort

No cable between the patient and monitor, so the patient is free from the bed. The transmitter is compact and light-weight so it is easy to wear.

## ❖ More efficient patient care

Instead of visit and manual measurement, automatic measurement saves the care staff's time.

## ❖ Easy to use

No complicated preparation to start measurement. The CALL button activates the recorder and nurse call at the central monitor. (ZS-910/920/930PK) The long battery life reduces running costs.

## Measurement Parameters



**ZS-910PK**  
ECG



**ZS-920PK**  
ECG, Resp



**ZS-930PK**  
ECG, Resp, SpO<sub>2</sub>

\*LCD shows SpO<sub>2</sub>



**ZS-940PK**  
NIBP, ECG, Resp, SpO<sub>2</sub>

\*LCD shows NIBP, SpO<sub>2</sub> and pulse rate

\*NIBP measurement mode can be selected from manual, continuous and interval

\*8-lead ECG measurement available(ZM-940PK)

## Options / Consumables

For full list of options/consumables, see the Technical Data separately available.

### Electrode lead



**K910A** 3-electrodes, snap

**K911** 3-electrodes, clip

### NIBP cuff



Reusable type with transmitter holder 13, 15 cm width

### SpO<sub>2</sub> probe (BluPRO®)



Reusable type

**P225H** finger, 0.6 m cord

**P225F** finger, 1.6 m cord

**P225G** multi-site



Disposable type Adult to child, 4 types

## Related Products



**Central monitor**  
CNS-9701J/K  
**Multiple patient receiver,** ORG-9700J/K



**Telemetry system**  
WEP-4204/4208J/K

This brochure may be revised or replaced by Nihon Kohden at any time without notice.

**NIHON KOHDEN**

**NIHON KOHDEN CORPORATION**

1-31-4 Nishiochiai, Shinjuku-ku, Tokyo 161-8560, Japan

Phone +81 (3) 5996-8036 Fax +81 (3) 5996-8100

www.nihonkohden.com